



# Developing an Assay to Determine the Signaling Pathway of GPCRs

Ashkhan Hojati, Takeharu Kawano, and Diomedes E. Logothetis




School of Medicine  
SURP Program  
Physiology and Biophysics Department  
Richmond, VA 23298

## Abstract

The G alpha 16 subunit is from the Gq protein family, but has been shown to be promiscuous (it can sometimes bind to Gi or Gs signaling receptors as well). Due to their promiscuity and intrinsic Gq behavior, they have been used to engineer chimeras that couple to specific receptors. By replacing varying lengths of their c-terminal ends with those of other G proteins' which couple to other pathways (e.g. Gi or Gs), chimeras can be constructed for further experimentation [1]. In this project the c-terminus 11-aa of human G alpha 16 was replaced with those of other proteins to construct G alpha 16 chimeras. These G alpha 16 chimeras can be particularly useful as they can couple to a receptor that couples to Gi or Gs for example while at the same time exhibiting Gq behavior (G alpha q signals PLC → hydrolysis of PIP2 to DAG and IP3 → IP3 mediated release of calcium from the endoplasmic reticulum → activates calcium-activated chloride channels). The activation of these chloride channels can be qualitatively seen by the “chloride spikes” which are ordinarily detected in Gq signaling using the Two-Electrode Voltage Clamp (TEVC). The hydrolysis of PIP2 from the membrane also results in channel current inhibition. This current inhibition can also be seen on channel recordings using the TEVC.

As previously mentioned, Gq signaling results in the increase of intracellular calcium after its release from the endoplasmic reticulum. This increase of intracellular calcium can also be measured in cells. Using calcium dyes, and a potentially working chimera we now have the tools to develop a high throughput assay in cells. This high throughput assay could be very beneficial as a system to screen receptors and drugs.

## Background



$G\alpha_{16}/\alpha i$   
 $G\alpha_{16}/\alpha s$

↑ [Ca]<sub>i</sub>

G protein-coupled receptors (GPCRs) have seven transmembrane helices, a C terminus located on the cytoplasmic side, and an N terminus on the extracellular side. GPCRs are involved in initiating many different cellular transduction mechanisms in our cells, and many present day drugs on the market (~50%) target GPCRs [3].

There is a critical interface between the receptor and the alpha subunit of the G protein. There are defined regions on G alpha subunit that interact with the receptor which couples to these proteins. Varying lengths of the c-terminal end of the alpha subunit are important in this interface.

